

Cytology Brushes

Single Use



Brush with Metal Tip

Efficient tube closure to secure smear yield

Set Versions Including Brush and Scissor

Safe and user-friendly specimen handling

Cytology Brushes – Single Use



Cytology Brush – Bronchoscopy

With conical metal tip and ergonomic handle

Catalogue No.	Description	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
GCB-01-18-120	Cytology Brush	120	1.8	2.0
GCB-01-24-120	Cytology Brush	120	2.3	2.8
GCB-02-18-120	Set (Brush & Scissor)	120	1.8	2.0



Cytology Brush – Gastroscopy

With conical metal tip and ergonomic handle

Catalogue No.	Description	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
GCB-01-18-180	Cytology Brush	180	1.8	2.0
GCB-01-24-180	Cytology Brush	180	2.3	2.8
GCB-02-18-180	Set (Brush & Scissor)	180	1.8	2.0



Cytology Brush – Colonoscopy

With conical metal tip and ergonomic handle

Catalogue No.	Description	Working Length (cm)	Sheath Outer Diameter (mm)	Min. Accessory Channel (mm)
GCB-01-18-230	Cytology Brush	230	1.8	2.0
GCB-01-24-230	Cytology Brush	230	2.3	2.8



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
 Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
 Email: sales@medi-globe.de • www.medi-globe.de



L-GCB-SU-00-EN / LL-0014 Vers. 8.0 (12/2020)

