

Stone Extraction Balloons

Single Use



L-GEB-SU-01-EN / LL-0033 Vers. 8.0 (12/2020)

Multi-Stage Extraction Balloons

Flexible diameter adjustment from 9–16 mm
in one balloon

Latex-Free Balloon Quality

Reduced risk of allergic reactions

Versions with Contrast Injection Above or Below the Balloon



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

Stone Extraction Balloons – Single Use




Single-Stage Stone Extraction Balloon – Double Lumen

Incl. syringe for max. air capacity; latex-free

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Contrast Injection	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-05-17-025	.025"	5	Above	12	1.5	200	2.8
GEB-05-23-035	.035"	7	Above	16	2.5	200	2.8


Stone Extraction Balloons – Single Use



Multi-Stage Stone Extraction Balloon - Triple Lumen

Incl. 3 syringes for max. air capacity; latex-free

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Contrast Injection	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-07-25-035	.035"	7.5	Above	9-13-16	1.5-2.5-5.0	200	2.8
GEB-07-25-035P	.035"	7.5	Below	9-13-16	1.5-2.5-5.0	200	2.8



Single-Stage Stone Extraction Balloon – Triple Lumen

Incl. syringe for max. air capacity; latex-free

Catalogue No.	Compatible with Guide Wire (inch)	Catheter Outer Diameter (French)	Contrast Injection	Balloon Diameter in Inflated Condition (mm)	Air Capacity (ml)	Working Length (cm)	Min. Accessory Channel (mm)
GEB-06-25-035	.035"	7.5	Above	16	5.0	200	2.8
GEB-06-25-035P	.035"	7.5	Below	16	5.0	200	2.8