

PolyCatch™ Polyp & Foreign Body Retrieval Devices

Single Use



Easy-opening Oval Shaped Retrieval Pouch

Efficient capture and retrieval of polyps
and foreign bodies

Rotatable Pouch

User friendly and precise device positioning

Latex-free Pouch

Avoids potential risk of latex allergic reactions

PolyCatch™ Retrieval Devices – Single Use



PolyCatch™ Polyp and Foreign Body Retrieval Device

Oval thermo-formed plastic bag

Catalogue No.	Bag Width (mm)	Bag Depth (mm)	Catheter Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GFB-08-15-120	15	12	1.8	120	2.0
GFB-08-25-180	25	17	2.6	180	2.8
GFB-08-25-230	25	17	2.6	230	2.8
GFB-08-35-230	35	20	2.6	230	2.8



L-GFB-SU-01-EN / LL-0111 Vers. 5.0 (01/2021)



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

