

Spray Catheters

Single Use



Application of Pigmented Liquids
for Marking of Suspicious Mucosa

Spray Catheters – Single Use



Wide Spray

Spray Catheter

Incl. stylet

Catalogue No.	Use for	Spray Head Diameter (mm)	Catheter Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSC-01-25-180	Wide Spray	2.4	2.4	180	2.8
GSC-01-25-230	Wide Spray	2.4	2.4	230	2.8



L-GSC-SUJ-00-EN / LL-0208 Vers. 2.0 (12/2020)



Medi-Globe GmbH

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392
Email: sales@medi-globe.de • www.medi-globe.de

