

# Injection Needles

Single Use



L-GSN-SU-01-EN / LL-0170 Vers. 2.0 (12/2020)

## **Eco-Injectra With Metal Tip Protection**

Reduces the risk of sheath puncture and secures a safe needle guidance in angulated endoscope positions

## **Ergonomic Handle & Needle Lock Mechanism**

Controlled puncture in targeted tissues and vessels

## **Transparent Needle Tube**

Flow control of therapeutics and fluids



### **Medi-Globe GmbH**

Medi-Globe-Straße 1–5 • D-83101 Rohrdorf OT Achenmühle  
Tel.: +49 (0) 8032 973-379 • Fax: +49 (0) 8032 973-392  
Email: [sales@medi-globe.de](mailto:sales@medi-globe.de) • [www.medi-globe.de](http://www.medi-globe.de)

# Injection Needles – Single Use



## Bronchoscopy - ECO-Injectra

Metal tip protection; needle locking mechanism (one hand operation)

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-21-18-120	0.5/25	4	1.8	120	2.0
GSN-22-18-120	0.5/25	6	1.8	120	2.0



## Gastroscopy - ECO-Injectra

Metal tip protection; needle locking mechanism (one hand operation)

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-21-23-180	0.5/25	4	2.3	180	2.8
GSN-22-23-180	0.5/25	6	2.3	180	2.8
GSN-23-23-180	0.7/22	4	2.3	180	2.8
GSN-24-23-180	0.7/22	6	2.3	180	2.8
GSN-25-23-180	1.0/19	4	2.3	180	2.8
GSN-26-23-180	1.0/19	6	2.3	180	2.8



# Injection Needles – Single Use



## Colonoscopy - ECO-Injectra

Metal tip protection; needle locking mechanism (one hand operation)

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-21-23-230	0.5/25	4	2.3	230	2.8
GSN-22-23-230	0.5/25	6	2.3	230	2.8
GSN-23-23-230	0.7/22	4	2.3	230	2.8
GSN-24-23-230	0.7/22	6	2.3	230	2.8
GSN-25-23-230	1.0/19	4	2.3	230	2.8
GSN-26-23-230	1.0/19	6	2.3	230	2.8



# Injection Needles – Single Use



## Bronchoscopy - Injectra

Removable safety clip; Luer-Lock-mechanism

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-04-23-120	1.0/19	6	2.3	120	2.8
GSN-03-19-120	0.7/22	6	1.9	120	2.2



## Gastroscopy - Injectra

Removable safety clip; Luer-Lock-mechanism

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-02-23-180	1.0/19	4	2.3	180	2.8
GSN-04-23-160	1.0/19	6	2.3	160	2.8
GSN-04-23-180	1.0/19	6	2.3	180	2.8
GSN-01-19-180	0.7/22	4	1.9	180	2.2
GSN-03-19-160	0.7/22	6	1.9	160	2.2
GSN-03-19-180	0.7/22	6	1.9	180	2.2
GSN-05-19-180	0.5/25	4	1.9	180	2.2
GSN-06-19-180	0.5/25	6	1.9	180	2.2



# Injection Needles – Single Use



## Colonoscopy - Injectra

Removable safety clip; Luer-Lock-mechanism

Catalogue No.	Needle Diameter (mm/Gauge)	Needle Length (mm)	Sheath Outer Diameter (mm)	Working Length (cm)	Min. Accessory Channel (mm)
GSN-02-23-230	1.0/19	4	2.3	230	2.8
GSN-04-23-230	1.0/19	6	2.3	230	2.8
GSN-01-19-230	0.7/22	4	1.9	230	2.2
GSN-03-19-210	0.7/22	6	1.9	210	2.2
GSN-03-19-230	0.7/22	6	1.9	230	2.2
GSN-05-19-230	0.5/25	4	1.9	230	2.2
GSN-06-19-230	0.5/25	6	1.9	230	2.2



